

At Laser-Laboratorium Göttingen, the performance of optical components for the UV and VUV spectral region is comprehensively characterized, using testing procedures in accordance with ISO standards:

*Absorptance (ISO 11551)*

*Scattering*

- total (ISO 13696)

- angular resolved

*Transmittance / Reflectance*

*Laser-induced damage threshold (ISO 11254)*

*Degradation (color centers)*

*Fluorescence / luminescence*

*Wavefront distortion (lens heating, compaction)*

## Employed radiation sources:

High power excimer lasers

**248nm, 193nm, 157nm**

Nd:YAG lasers

**1064nm, 355nm, 266nm**

laser-induced EUV source

**13nm**

## Absorptance

Absolute absorptance is determined using a high-resolution laser-calorimetric technique with sensitivity in the ppm range. Both linear and non-linear (two-photon) absorption losses can be recorded. Measurements permit a rapid assessment of the long-term performance of DUV/VUV optics at low energy density levels.

*direct measurement of absolute absorptance in accordance with ISO 11551*

*determination of linear and non-linear absorption coefficients*

*degradation rates without thermal loading of the samples*

*non-destructive quality assessment*

## Fluorescence

Temporally-resolved laser-induced fluorescence / luminescence can be monitored during the absorptance measurements with the help of a gateable image-intensified UV spectrometer (cf. Fig. 1).




Fig. 1: Luminescence spectra of two different  $\text{CaF}_2$  samples irradiated at 157nm

## Scattering

Scatter losses in UV optics are determined using a Total Scattering (TS) setup with an excimer laser as pulsed light source and a high-precision Coblentz hemisphere as integrating element, allowing to monitor both forward and backward scattering independently. Moreover, angular resolved scatter data can be obtained.

*'at wavelength' recording of scatter losses in the DUV spectral range*

*sensitivity: < 1 ppm at 248nm*

*< 30 ppm at 193nm*

*scatter maps (cf. Fig. 2, lateral resolution 0.5mm)*




Fig. 2: Scatter map of super-polished  $\text{CaF}_2$  at 248nm

## Wavefront monitoring during UV irradiation:

Permanent or transient laser-induced deformation of DUV/VUV optical elements (e.g. due to lens-heating, compaction or rarification in fused silica) can be monitored on-line with a resolution in the nm range using a Hartmann-Shack wavefront sensor.




Fig. 3: Compaction of fused silica during ArF laser irradiation monitored with Hartmann-Shack wavefront sensor; above: experimental set-up; left: before irradiation; right: during and after irradiation (densification  $\sim 10^{-6}$ ).